

CertoClav Connect

Fully automatic, with a 7-inch CertoTouch screen that explains every step on video.

Fully automatic laboratory autoclave with 7 inch high-res touch screen and intuitive user interface. All user tasks are explained with videos on the machine.



Technical data

Chamber volume	18 L
Chamber diameter	26 cm
Chamber depth	35 cm
Chamber material	Stainless steel 316L
Max. operating pressure	2.7 bar
Temperature	20 – 140 °C, adjustable
Heating power	1300 W
Electrical supply	230 V, 50/60 Hz
Net weight	16 kg
Dimensions (W × D × H)	44 × 38 × 54 cm
Orientation	Vertical, top loading
Control	Fully automatic
Fractionated vacuum	No
Flexible media probe	Yes
Interfaces	WiFi, LAN, USB, CertoCloud
Country of origin	Austria

Released for

STERILISATION OF

- Unwrapped solids without hollow bodies
- Liquids
- Bottles
- Instruments
- Plastics
- Rubber (105 °C)

PREPARATION OF

- Culture media and agar

WASTE DISPOSAL

- Sterilisation of waste
- Waste in disposal bags

TESTING AND CYCLES

- Material testing
- Long cycles, up to 7 days

FUNCTION

- Digital documentation
- Programmable
- Touchscreen
- CertoCloud

EQUIPMENT AND CONNECTION

- Standalone, no plumbing needed
- Single-phase socket, no three-phase supply

Loading capacity: Stainless steel series 18 L

Chamber Ø 260 mm, height 330 mm, effective height 290 mm.

LOAD CARRIER	ARTICLE	H × Ø MM	QTY
Instrument tray	1300155	41 × 250	1
Wire basket	1300149	250 × 250	1

BOTTLE	ML	H × Ø MM	BOTTLES
DURAN® Original GL 45	2000	265 × 136	1
DURAN® GLS 80	2000	248 × 136	1
DURAN® Original GL 45	1000	230 × 101	4
DURAN® GLS 80	1000	218 × 101	4
DURAN® Original GL 45	500	181 × 84	5
DURAN® Youtility	500	193 × 78	5
DURAN® GLS 80	500	148 × 101	4
DURAN® Original GL 45	250	143 × 70	9
DURAN® Youtility	250	158 × 66	9
DURAN® GLS 80	250	105 × 95	4
DURAN® Youtility	125	124 × 55	12
DURAN® Original GL 45	100	105 × 56	12

Source: CertoClav loading capacity sheet, 2022-02. DURAN® is a trademark of DWK Life Sciences.

Accessories

ACCESSORY	ARTICLE NUMBER
Instrument plate with feet	1252203
Service set 2	1310186
Wire basket, Ø 25 cm, with handle	1300149
Media probe	1300185
CertoCloud Premium licence	8500212
Water distiller	8500482

CE marked. Conforms to 2014/35/EU (Low Voltage Directive), 2014/30/EU (EMC Directive), 2011/65/EU (RoHS Directive). Figures apply to article 8510176; other versions may differ. Values marked "on request" are not published in this document.