




Beladungskapazität Loading Capacity CertoClav Aluminium-Serie 12 L CertoClav Aluminum Series 12L



**Kammer
Chamber**

**Instrumentenplatte
Instrument tray**
Artikelnr.: 1252203

**Drahtkorb
Wire basket**
Artikelnr.: 8583010



Artikel	Artikelnr.	Höhe/ Hight [mm]	Durchmesser/ Diameter [mm]	Volume
CertoClav EL 115/121 °C	8501312A	260	240	12 l
CertoClav EL 125/140 °C	8501312	260	240	12 l
CertoClav MultiControl	8502312A	260	240	12 l

Artikel	Artikelnr.	Höhe/ Hight [mm]	Durchmesser/ Diameter [mm]	
Instrument tray	1252203	41	225	1 Stk.
Wire basket	8583010	170	230	1 Stk.

Flasche / Bottle	Volumen / Capacity [ml]	Höhe / Height [mm]	Durchmesser / Diameter [mm]	Beladungskapazität / Loading capacity
DURAN® Original GL 45	1000	230	101	1 (3*)
DURAN® GLS 80	1000	218	101	1 (3*)
DURAN® Original GL 45	500	181	84	5
DURAN® Youtility	500	193	78	5
DURAN® GLS 80	500	148	101	3
DURAN® Original GL 45	250	143	70	7
DURAN® Youtility	250	158	66	8
DURAN® GLS 80	250	105	95	4
DURAN® Youtility	125	124	55	12
DURAN® Original GL 45	100	105	56	12

*nur/ only MultiControl



Beladungskapazität

Loading Capacity


CertoClav Aluminium-Serie 18L

CertoClav Aluminum Series 18L

Kammer
Chamber

Instrumentenplatte
Instrument tray
Artikelnr.: 1252203

Drahtkorb
Wire basket
Artikelnr.: 8583010



Artikel	Artikelnr.	Höhe/ Hight [mm]	Durchmesser/ Diameter [mm]	Volume
CertoClav EL 115/121 °C	8501318A	390	240	18 l
CertoClav EL 125/140 °C	8501318	390	240	18 l

Artikel	Artikelnr.	Höhe/ Hight [mm]	Durchmesser/ Diameter [mm]	
Instrument tray	1252203	41	225	1 Stk.
Wire basket	8583010	170	230	1 Stk.

Flasche / Bottle	Volumen / Capacity [ml]	Höhe / Height [mm]	Durchmesser / Diameter [mm]	Beladungskapazität / Loading capacity
DURAN® Original GL 45	2000	265	136	1
DURAN® GLS 80	2000	248	136	1
DURAN® Original GL 45	1000	230	101	3
DURAN® Youtility	1000	253	93	4
DURAN® GLS 80	1000	218	101	3
DURAN® Original GL 45	500	181	84	5
DURAN® Youtility	500	193	78	5
DURAN® GLS 80	500	148	101	3
DURAN® Original GL 45	250	143	70	7
DURAN® Youtility	250	158	66	8
DURAN® GLS 80	250	105	95	4
DURAN® Youtility	125	124	55	12
DURAN® Original GL 45	100	105	56	12