



Beladungskapazität

Loading Capacity


CertoClav Edelstahl-Serie

CertoClav Stainless Steel Series

Kammer
Chamber

Instrumentenplatte
Instrument tray
Artikelnr.: 1300155

Drahtkorb
Wire basket
Artikelnr.: 1300149



Artikel	Artikelnr.	Höhe/ Hight [mm]	Durchmesser/ Diameter [mm]	Volume
CertoClav Classic 18 l 115/121 °C	8510184	330	260	18 l
CertoClav Classic 18 l 125/140 °C	8510174	330	260	18 l
CertoClav MultiControl 2	8510175	330	260	18 l
CertoClav Connect	8510176	330	260	18 l

Artikel	Artikelnr.	Höhe/ Hight [mm]	Durchmesser/ Diameter [mm]	
Instrument tray	1300155	41	250	1 Stk.
Wire basket	1300149	250	250	1 Stk.

Flasche / Bottle	Volumen / Capacity [ml]	Höhe / Height [mm]	Durchmesser / Diameter [mm]	Beladungskapazität / Loading capacity
DURAN® Original GL 45	2000	265	136	1
DURAN® GLS 80	2000	248	136	1
DURAN® Original GL 45	1000	230	101	4
DURAN® GLS 80	1000	218	101	4
DURAN® Original GL 45	500	181	84	5
DURAN® Youtility	500	193	78	5
DURAN® GLS 80	500	148	101	4
DURAN® Original GL 45	250	143	70	9
DURAN® Youtility	250	158	66	9
DURAN® GLS 80	250	105	95	4
DURAN® Youtility	125	124	55	12
DURAN® Original GL 45	100	105	56	12